

















MYO/WIRE II™ Monofilament Sternotomy Wires

Needle (Shown Actual Size)		USP SIZE	MYO/Wire II™ Monofilament Sternotomy Wires (12 packages per box)			
Point Configuration	Needle Length		2 x 18" (45cm) Per Package Cat. No.	4 x 18" (45cm) Per Package Cat. No.	6 x 18" (45cm) Per Package Cat. No.	2 x 30" (75cm) Per Package Cat. No.
 <p>A&E Standard Taper Cutting</p>	 55mm 3/8 Circle	5	----	045-002	045-003	045-101*
		6	----	046-002	046-003*	046-101*
		7	----	047-002*	047-003	----
 <p>Conventional Cutting CCS#1</p>	 49mm 1/2 Circle	5	045-031*	045-032	045-033	045-131
		6	046-031	046-032	046-033	046-131
		7	047-031*	047-032	047-033	----
 <p>Conventional Cutting CCS#2</p>	 59mm 1/2 Circle	5	045-061*	045-062*	045-063*	045-161
		6	046-061	046-062	046-063*	046-161
		7	047-061	047-062	047-063	----
 <p>Conventional Cutting CCS#3</p>	 49mm 1/2 Circle	5	045-081*	045-082	045-083*	045-181*
		6	046-081*	046-082	046-083*	046-181*
		7	047-081*	047-082	047-083*	----
 <p>Diamond Point DP</p>	 48mm 1/2 Circle	5	----	045-042*	045-043*	045-141*
		6	----	046-042*	046-043*	046-141*
		7	----	047-042	047-043*	----
 <p>Reverse Cutting RCE#1</p>	 55mm 1/2 Circle	5	045-051*	045-052*	045-053*	045-151*
		6	046-051*	046-052*	046-053	046-151*
		7	047-051*	047-052*	047-053*	----
 <p>Taper Cutting TC#1</p>	 51mm 1/2 Circle	5	----	045-072	045-073	045-171*
		6	----	046-072	046-073	046-171*
		7	----	047-072	047-073	----
 <p>Taper Cutting TC#2</p>	 40mm 1/2 Circle	1	----	041-082	----	----
		2	----	042-082	----	----

*Non-stock product. Available as custom order product in minimum 10 box quantities.